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
The Foot Alignment Socks use instructions

The Foot Alignment Socks give comfort and help to relieve common foot ailments such as bunions, hammertoes, plantar fasciitis, foot cramps, toe cramps, crooked toes and general foot discomfort. They are mainly used as a **supplement to active exercises and physiotherapy**. Foot Alignment Socks improve blood circulation and give your feet and toes new energy. They help by keeping your feet and toes mobile, flexible and less painful. To really relieve your feet from pain, you should follow a few rules when wearing Socks.

Proceed slowly. If irritation occurs, discontinue use.

- 1) The embroidery feet should point to each other when the Socks are put on. Start with the dividers at the tip of your toes.
(wear for approx. 15-20 minutes per day, kids 5-10 minutes)
- 2) After a few days of wearing your Socks, slide the dividers half way down between your toes.
(wear for approx. 1 to 2 hours per day, kids 30-60 minutes)
- 3) Within a couple of weeks you should be able to slide the toe dividers all the way down between your toes providing maximum benefit and comfort to your feet (it is possible to wear the Socks overnight while sleeping now).
Sleeping with the Socks on is not recommended for children.

Foot Alignment Socks are NOT intended for shoes, nor for normal walking. They are **most effective when you are relaxing**, when the feet are not under pressure – i.e. when sitting, lying down or sleeping. **If the sleeping with Socks is uncomfortable for you, it is not required to achieve the benefits.** Follow your own feelings, take off the Socks if needed and continue wearing them the next day.

Fiber: 90 % cotton, 7 % polyester, 3 % elastane

Washing instructions: by hand or machine wash cold (30 °C), do not tumble dry

These socks are not a substitute for medical treatment. Solve foot problems in time with a physiotherapist.

Exclusive importer for Europe:

pronožky.cz

www.pro-nozky.cz/en, info@pro-nozky.cz

WEAR THE SOCKS CORRECTLY TO INCREASE THEIR EFFECT.

Adjustační ponožky®, Adjustky®, pro nožky.cz®, Adjustačné ponožky® and The Original Foot Alignment Socks® are registered trademarks with application numbers [497166](#), [581608](#), [516612](#), [763-2019](#) and [17979784](#). The shape of the socks is also protected by trademarks (no. [546682](#), [17979786](#) and [772-2019](#)). Original U.S.patent 7784115B1.